

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | AMG888, U3-1287 |
| Conjugate | PE-conjugated |
| Synonyms | HER3;ERBB3 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:400 |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Homo sapiens |
| IgG type | Human IgG1 - kappa |
| Reactivity | Human |
| Uniprot ID | P21860 |
| Target | HER3 |
| Description | PE-conjugated Anti-HER3 (patritumab biosimilar)mAb |
| Delivery | In Stock |
| Storage&Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |






Figure 1. Flow cytometry analysis with 100 μ l/test (1:400) PE-conjugated Anti-HER3 (patritumab biosimilar)mAb (BME100057P) on HEK293 cells transfected with human HER3 (Blue histogram) or HEK293 cells (Red histogram).

DIMA BIO CONFIDENTIAL

