

PRODUCT INFORMATION

|                              |                                                                                                                                                                                                                                                                        |
|------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clone ID                     | 4A12                                                                                                                                                                                                                                                                   |
| Target                       | Human IgG                                                                                                                                                                                                                                                              |
| Synonyms                     | N/A                                                                                                                                                                                                                                                                    |
| Host Species                 | Rabbit                                                                                                                                                                                                                                                                 |
| Description                  | PE-conjugated Anti-Human IgG antibody(4A12), Rabbit mAb                                                                                                                                                                                                                |
| Delivery                     | Under Development                                                                                                                                                                                                                                                      |
| Uniprot ID                   | N/A                                                                                                                                                                                                                                                                    |
| IgG type                     | Rabbit IgG                                                                                                                                                                                                                                                             |
| Clonality                    | Monoclonal                                                                                                                                                                                                                                                             |
| Reactivity                   | Human                                                                                                                                                                                                                                                                  |
| Applications                 | Flow Cyt                                                                                                                                                                                                                                                               |
| Recommended Dilutions        | Flow Cyt 1:100                                                                                                                                                                                                                                                         |
| Purification                 | Purified from cell culture supernatant by affinity chromatography                                                                                                                                                                                                      |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA                                                                                                                                                                                                                               |
| Storage & Shipping           | Store at 2°C-8°C for 6 months                                                                                                                                                                                                                                          |
| Background                   | N/A                                                                                                                                                                                                                                                                    |
| Usage                        | Research use only                                                                                                                                                                                                                                                      |
| Conjugate                    | PE-conjugated                                                                                                                                                                                                                                                          |
| DIMA Disclaimer              | All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. |

