

PRODUCT INFORMATION

Common Name	RO4858696-000,teprotumumab-trbw
Conjugate	PE-conjugated
Synonyms	CD221, IGFIR, IGFR, JTK13
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	Human IgG1(E356D,M358L) - kappa
Reactivity	Human
Target	IGF-1R
Uniprot ID	P08069
Description	PE-conjugated Anti-IGF-1R(teprotumumab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

