

PRODUCT INFORMATION

Common Name	REGN668, SAR23189
Conjugate	PE-conjugated
Synonyms	CD124;IL-4RA;IL4RA
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	Human IgG4 - Kappa
Reactivity	Human
Target	IL4RA
Uniprot ID	P24394
Description	PE-conjugated Anti-IL-4RA (dupilumab biosimilar)mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

