

PRODUCT INFORMATION

Common Name	AMG827, KHK4827
Conjugate	PE-conjugated
Synonyms	CD217;CDw217;IL-17RA;IL17R;CANDF5;hIL-17R
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	Human IgG2 - Kappa
Reactivity	Human
Target	IL17RA
Uniprot ID	Q96F46
Description	PE-conjugated Anti-IL17RA (brodalumab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

