

PRODUCT INFORMATION

Common Name	MABp1
Conjugate	PE-conjugated
Synonyms	IL-1 alpha;IL-1A;IL1;IL1-ALPHA;IL1F1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	IL1A
Uniprot ID	P01583
Description	PE-conjugated Anti-IL1A(bermekimab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

