

PRODUCT INFORMATION

Common Name	ANB-019
Conjugate	PE-conjugated
Synonyms	ILRL2; IL-36R; IL1RRP2; IL-36R
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric Humanized
IgG type	Human IgG4 - Kappa
Reactivity	Human
Target	IL1RL2
Uniprot ID	Q9HB29
Description	PE-conjugated Anti-IL1RL2(imsidolimab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

