

PRODUCT INFORMATION

Uniprot ID	P01589
Common Name	CHI621
Conjugate	PE-conjugated
Synonyms	IL2RA;CD25;p55;IL2-RA;IL-2-RA
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	IL2RA
Description	PE-conjugated Anti-IL2RA (basiliximab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

