

PRODUCT INFORMATION

Uniprot ID	P16871
Common Name	Effi-7, Effi7-h, MD-707
Conjugate	PE-conjugated
Synonyms	CD127; CDW127; IL-7R-alpha; IL7RA; ILRA
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG4
Reactivity	Human
Target	IL7RA
Description	PE-conjugated Anti-IL7RA(lusvertikimab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

