

PRODUCT INFORMATION

Common Name	7F3com-2H2, MEDI-528
Conjugate	PE-conjugated
Synonyms	Interleukin-9;IL-9;Cytokine P40;T-cell growth factor P40
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	IL9
Uniprot ID	P15248
Description	PE-conjugated Anti-IL9(enokizumab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

