

PRODUCT INFORMATION

Common Name	AN100226,AN100226,BG00002
Conjugate	PE-conjugated
Synonyms	CD49D;IA4
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG4 - Kappa
Reactivity	Human
Target	ITGA4
Uniprot ID	P13612
Description	PE-conjugated Anti-ITGA4 (natalizumab biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

