

**PRODUCT INFORMATION**

<b>Common Name</b>	TBC
<b>Conjugate</b>	PE-conjugated
<b>Synonyms</b>	hK2;hGK-1;KLK2A2
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	1µl/test
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with,0.18%BSA□0.1%Proclin300
<b>Host Species</b>	N/A
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	KLK2
<b>Uniprot ID</b>	P20151
<b>Description</b>	PE-conjugated Anti-KLK2(alpitatug biosimilar) mAb
<b>Delivery</b>	In Stock
<b>Storage&amp;Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

DIMABIO CONFIDENTIAL

