

PRODUCT INFORMATION

Common Name	MK-0482-Hu52B8
Conjugate	PE-conjugated
Synonyms	ILT3; LIR5; CD85k
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid-PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG4 - Kappa
Reactivity	Human
Target	LILRB4
Uniprot ID	Q8NHJ6
Description	PE-conjugated Anti-LILRB4(MK-0482 biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

