

PRODUCT INFORMATION

| | |
|------------------------------|--|
| Uniprot ID | P09603 |
| Common Name | MCS-110, MCS110 |
| Conjugate | PE-conjugated |
| Synonyms | CSF-1;MCSF |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Humanized |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | M-CSF |
| Description | PE-conjugated Anti-M-CSF(lacnotuzumab biosimilar) mAb |
| Delivery | Under Development |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| DIMA Disclaimer | All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scr |

