

PRODUCT INFORMATION

Uniprot ID	P15941
Common Name	PAB001-PAB001-ADC-hMUC1-1H7
Conjugate	PE-conjugated
Synonyms	MUC-1; CD227; PUM
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid-PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1(E239D, M241L)
Reactivity	Human
Target	MUC1
Description	PE-conjugated Anti-MUC1(PAb-001 biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

