

PRODUCT INFORMATION

Uniprot ID	P26718
Common Name	5C5
Conjugate	PE-conjugated
Synonyms	NKG2D;CD314;KLRK1;NK cell receptor D
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	NKG2D
Description	PE-conjugated Anti-NKG2D(5C5 biosimilar, L234A/L235A) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

