

PRODUCT INFORMATION

Uniprot ID	Q13308
Common Name	PF-06523435,hu24
Conjugate	PE-conjugated
Synonyms	CCK4
Applications	Flow Cyt
Recommended Dilutions	1μl/test
Formulation & Reconstitution	Liquid□PBS with,0.18%BSA□0.1%Proclin300
Host Species	Humanized
IgG type	IgG1(E239D, M241L)
Reactivity	Human
Target	PTK7
Description	PE-conjugated Anti-PTK7(cofetuzumab biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 1µl/test PE-conjugated Anti-PTK7(cofetuzumab biosimilar) mAb on 293T cells (Blue histogram) or Raji cells (Red histogram).

