

PRODUCT INFORMATION

Common Name	JR-141
Conjugate	PE-conjugated
Synonyms	TR;TfR;TfR1;Trfr;T9;p90;CD71
Applications	Flow Cyt
Endotoxin	Less than 1.0 EU/ μ g by the LAL method. For <1 EU/mg requirements, please contact us for customization.
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	TFRC
Uniprot ID	P02786
Description	PE-conjugated Anti-TFRC(pabinafusp alfa biosimilar) mAb
Delivery	Under Development
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

