

PRODUCT INFORMATION

Common Name	AGX101
Conjugate	PE-conjugated
Synonyms	M3S1;TAAL6
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1µl/test
Formulation & Reconstitution	Liquid□PBS with,0.18%BSA□0.1%Proclin300
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	TM4SF1
Uniprot ID	P30408
Description	PE-conjugated Anti-TM4SF1 (AGX101 biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 1 μ l/test PE-conjugated Anti-TM4SF1 (AGX101 biosimilar) mAb on SNU-5 cell line (Blue histogram) or 293T cell line (Red histogram).

DIMA BIO CONFIDENTIAL

